

Video acquisition with JAI Cameras

The **videoJaiAq** toolbox has been developed by F. Ciriello to provide a framework for G. R. Hunt's research group to share acquisition and processing tools. It uses the **videoJaiAq** class, which you can download from the **accumulator** (or on request at **fc397@cam.ac.uk**) and added to your MATLAB userpath. You can see your user path by running

```
userpath
```

in the command prompt.

Make sure you have updated MATLAB! The toolbox has been tested on MATLAB 2019a.

I would recommend to set the **current directory** to the main project folder you are working in to ensure that outputs are stored in an appropriate location.

```
cd D:\Users\fc397\Documents\PROJECT_Flow_in_a_corner\myProject
```

Contents

How to use the JAI video acquisition toolbox.....	1
Acquisition functionalities.....	2
Basics.....	2
Load camera.....	2
Camera controls.....	2
Preview.....	3
Acquisition.....	4
Playback.....	4
Preview Tools.....	4
Cropping.....	4
Calibration.....	5
Multi-camera acquisition.....	6
Stereo calibration.....	6
Live image processing.....	6
Background Subtraction.....	6
Post-processing.....	7
Installing the JAI cameras to use in MATLAB.....	7
Known compatibility issues.....	8

How to use the JAI video acquisition toolbox

I have created the video acquisition class (**videoJaiAq**) to process video streams from monochrome scientific cameras that use the *GenICam* protocol (tested on JAI Spark 5000-M-USB).

Users can use **help** functions for more information on any of the functions used in this class.

For example, to get more info on the **videoJaiAq** class (as a whole) by simply typing in the command window

```
help videoJaiAq
```

The class allows a user to create an "*acquisition*" object and use in-built methods to perform common tasks, including:

- load a camera videostream (**videoJaiAq**)
- preview a live stream from the camera (**preview**, **advancedPreview**)
- crop to a region of interest ROI (**crop**)
- open a simplified GUI to simplify for camera settings (**camBasicSettings**)
- change, save and load pre-defined settings (**camSettings**)
- recolor live stream in false colour (**camColor**)
- open a tool to assist in manual focus (**camFocus**)
- calibrate camera (**camCalibrate**)
- acquire video and save into a .avi file (**acquire**)
- acquire background frames and save into a .avi movie and .tif image (**acquireBackground**)

Furthermore, the toolbox can be used to apply a range of **live image processing** algorithms, including:

- background subtractions (**livebackgroundSubtract**)
- detect edges (**liveEdge**)
- motion detection with optical flow methods (**liveOF**)
- process singleshot PIV image pair (**livePIV**)

and post-processing features:

- extract raw frames from movie file (**extractFrames**)
- background subtract images (**backgroundSubtract**)

Acquisition functionalities

In case of trouble, reset camera connection via:

```
imaqreset;
```

Basics

Load camera

Create an acquisition object (**cam1**) using the **videoJaiAq** class file. The first argument allows you to choose which camera to connect to.

```
cam1 = videoJaiAq(1, 'Mono8');
```

Most functionalities have been tested on 8-bit depth, but I will work on extending to 10 and 12 bit depths.

Camera controls

The **camBasicSettings** function opens up a dialog box with the most commonly-used options for camera settings.

```
camBasicSettings(cam1)
```

Note that **camBasicSettings** cannot be used when **preview / acquisition** is active if *exposure fps = acquisition fps*. Also note that **camBasicSettings** does not work with the **advancedPreview** function due to callback errors.

All other controls can be accessed simply via **cam1.camera.src.controlName**. For example, to change the gain of the camera you can use:

```
cam1.camera.src.Gain = 2;
```

A list of camera controls can be inspected with the **get** function.

```
get(cam1.camera.src)
```

The **camSettings** function adds save and load functionalities. It also has a wider selection of options for camera settings.

```
camSettings(cam1)
camSettings(cam1, 'save', 'mySettings')
camSettings(cam1, 'load', 'mySettings')
```

Preview

The live stream from the camera is contained with the **obj.liveStream** property handle. You can modify the live stream in different ways. Axis properties are in **obj.liveStream.Parent** and figure properties are in **obj.liveStream.Parent.Parent**.

Note that the **preview** function from the image acquisition toolbox has been overloaded for that it can take in the **videoJaiAq** class and it creates its own figure handle object (*i.e.* the aforementioned **obj.liveStream**) from the videoinput object.

```
preview(cam1)
```

An **advancedPreview** function is also available to inspect the grey levels in the image and detect saturation levels.

```
advancedPreview(cam1);
```

The **camColor** function allows to change the colormap for the current live stream. The live stream must be open. This is a useful feature to detect saturation levels. I have created a bespoke colour map to highlight saturated regions called **satColor** (a grey scale with red and green indicating black and white saturated respectively).

```
camColor(cam1, gray)    % default
camColor(cam1, jet)     % you can use any pre-defined MATLAB colormap
camColor(cam1, 'satColor')
```

Acquisition

The **acquireBackground** function acquires 100 frames by default and stores the average of the acquired frames as the background in **cam1.background.image**. An optional second argument can be used to change the number of background frames to acquire.

```
acquireBackground(cam1,10)
```

You can inspect the background image with:

```
imshow(cam1.background.frame)
```

The **acquire** function can be used to record a avi movie. A dialog box appears to start acquisition manually.

```
acquire(cam1,'on','on',10,'myMovieName')
```

Default outputs are named '**raw.avi**' and placed in the **exp*** folder.

A second optional argument **saves the background video** to the same directory ('on' or 'off').

A third optional argument **saves the calibration** to the same directory ('on' or 'off').

A fourth optional argument specifies the **number of frames** to acquire for.

A fifth optional argument manually sets the **filename**.

A sixth optional argument specifies whether to **acquire with live stream** 'on' or 'off'.

A single image can be acquired with

```
img = getsnapshot(cam1.camera.vid);  
imwrite(img, 'outputs/image.tif');
```

Playback

Any video can be replayed using MATLAB's **implay** function.

```
implay('outputs/exp1/raw.avi')
```

Preview Tools

Cropping

The **crop** function can be used after **preview** or as a standalone **preview&crop** function.

The **advancedPreview** function can also be used in conjunction with **crop**. The function asks for a user-input for an ROI that can be drawn directly on the live stream. The ROI can be adjusted and confirmed with a double-click.

```
crop(cam1)
```

The **cropReset** function resets the cropping to the original camera resolution.

```
cropReset(cam1)
```

The **camFocus** function has been designed to assist in manual focus mode. A user-prompted rectangle is required to zoom into a region that requires focus. This opens up a new figure window with the zoomed in image and a plot that shows strips of pixel intensities. The more detail in the single the better the focus.

```
camFocus(cam1);
```

Calibration

The **camCalibrate** function assists calibration of the camera by taking processing images of a real-world object of known size. All properties stored in **obj.calibration**.

A '**line**' argument applies a simple linear calibration of px to mm requiring a user-input.

```
camCalibrate(cam1,'line');
```

A '**checker**' board argument guides the user through a more advanced calibration. x images (by default, $x = 10$, at least 3, ideally $10 - 20$) are acquired of a checkerboard placed at different orientations and positions. The function deduces *extrinsic* camera parameters (**obj.calibration.cameraParams**). To simplify use, a px size in mm is stored in **obj.calibration.pxsize** relative to the first calibration plane.

```
camCalibrate(cam1,'checker',squareSize,numCalibrationImages) % squareSize in mm  
camCalibrate(cam1,'checker',13,5)
```

The checker board calibration, amongst other functionalities, allows you to reproject images so that you can remove warping effects using the **undistortImage** function:

```
rawImage = getsnapshot(cam1.camera.vid);  
undistortedImage = undistortImage(rawImage,cam1.calibration.cameraParams);  
subplot(1,2,1)  
imshow(rawImage)  
subplot(1,2,2)  
imshow(undistortedImage)
```

For GUI guided camera calibration, you can use MATLAB's **cameraCalibrator** app:

```
cameraCalibrator
```

and find **more information** on camera calibration: <https://uk.mathworks.com/help/vision/ug/camera-calibration.html>

A '**dye**' argument can be used to apply a dye calibration using a dye-filled prism.

```
camCalibrate(cam1,'dye') % this feature is under development
```

Multi-camera acquisition

Several cameras can be loaded as different objects

```
cam1 = videoJaiAq(1, 'Mono8');  
cam2 = videoJaiAq(2, 'Mono8');
```

Multi-camera functionalities are the same as the single-camera functionalities described in the **Basics** section. Simply insert the camera objects as an array.

```
preview([cam1, cam2])  
acquireBackground([cam1, cam2],10)  
acquire([cam1, cam2], 'on', 'on', 100, 'myMovieName', 'on');
```

Stereo calibration

The **camStereoCalibrate** function can be used to calibrate to or more cameras to enable stereoscopic vision functionalities. ♦

```
camStereoCalibrate([cam1, cam2],13,10) % args = {[camarray],squareSize,numCalibImages}
```

MATLAB's GUI for stereo calibration can be launched with:

```
stereoCameraCalibrator
```

Once the calibration is complete, you can acquire stereoscopic images and videos.

```
rawImage1 = getsnapshot(cam1.camera.vid);  
rawImage2 = getsnapshot(cam2.camera.vid);  
  
[rectImage1, rectImage2] = rectifyStereoImages(rawImage1, rawImage2, cam1.stereoCalibration.st  
  
disparityMap = disparity(rectImage1, rectImage2);  
  
imshow(stereoAnaglyph(rectImage1, rectImage2))  
imagesc(disparityMap)
```

Live image processing

There are a number of functionalities here that I am developing.

Background Subtraction

A background subtraction can be applied to the live stream by using:

```
liveBackgroundSubtract(cam1)
```

Post-processing

The `extractFrames` function can be used to batch save all images from a file.

```
extractFrames('outputs/exp1/raw.avi')
```

I would recommend not to use **extractFrames**. You can read frames using a **videoreader** object:

```
video = VideoReader('outputs/exp1/raw.avi')  
frame = read(video,2)    % reading second frame (n.b. first frame often black)
```

The **makeVideo** function opens a GUI with **imaqmontage** and allows you to choose background subtraction, fps, moviename, and contrast adjustments.

Installing the JAI cameras to use in MATLAB

There are a few necessary steps to setup acquisition with the JAI camera in MATLAB and use the toolbox. These are:

1. Install **Jai SDK** from <https://www.jai.com/support-software/jai-software>.
2. Install MATLAB with the **image acquisition**, **image processing**, **computer vision** toolboxes from <https://uk.mathworks.com/downloads/>.
3. Install MATLAB hardware **adaptor for GenICam** (search for 'genicam adaptor' in MATLAB Apps)
4. Add Jai SDK directory to **environment variables**

go to **System Properties > Advanced > Environment Variables**

add new system variable

> Variable Name : **GENICAM_GENTL64_PATH**

> Variable value : **C:\Program Files\JAI\SDK\bin**

Change the variable name depending on whether computer is 32 or 64 bits and the variable value depending on your path directory for the Jai SDK. I would suggest to either type the variable names and values or to copy & paste these in after you got rid of formatting (e.g. you could do this by copy & pasting these in your browser search bar and then copy & paste again in the environment variable dialog boxes).

5. Restart MATLAB

Known compatibility issues

1. **False colour:** Older versions of MATLAB (prior 2017b) hold the **colormap** property in the figure object rather than the axis. This causes problems when using false color in live streams.