Methods	Description
Wethous	Description
<pre>cam = videoJaiAq(1, 'Mono8');</pre>	load object called cam
camBasicSettings(cam);	dialog box – simple settings
camSettings(cam);	dialog box – extended settings
<pre>camSettings(cam, 'save', 'filename');</pre>	save settings to cameraSettings folder
<pre>camSettings(cam, load, 'filename');</pre>	load settings from cameraSettings folder
<pre>cam.camera.src.PropertyName = PropertyValue;</pre>	change camera property
get(cam.camera.src)	list camera properties
preview(cam)	preview camera
advancedPreview(cam)	preview camera with live histogram
crop(cam)	crop to user-defined ROI
cropReset (cam)	reset crop
camFocus (cam)	open focus assist tool
Camp Ocus (Cam)	open rocus assist tool
camColor(cam,gray)	change colour of live stream
<pre>camColor(cam, jet) camColor(cam, 'satColor')</pre>	
Camcolor (Cam, Sattorn)	
acquireBackground(cam, x)	acquire x number of background images
	and save a time-average
<pre>camCalibrate(cam, 'line') camCalibrate(cam,'checker', squareSize, numCalImages) camCalibrate(cam, 'dye')</pre>	open user-assisted calibration
acquire(cam, 'on', 'on', x,'myMovieName', 'on')	acquire x images and save as avi movie in outputs/exp* folder
	2 nd argument sets background save status
	3 rd argument sets calibration save status
	6 th argument sets live acquisition status
<pre>fieldnames(cam.background) fieldnames(cam.calibration)</pre>	display field names of a structure
<pre>i = getsnapshot(cam.camera.vid);</pre>	and simple income
imshow(i)	get single image
<pre>imwrite(i, 'filename')</pre>	show image save image
<pre>imread('filename')</pre>	load image
	Todu IIIage
<pre>implay('outputs/exp*/raw.avi')</pre>	play movie
liveBackgroundSubtract(cam)	live background subtraction mode
liveEdge (cam)	live edge detection
livePIV(cam)	process a PIV image pair
liveOF(cam)	live optical flow mode
TIVOOL (Com)	nve optical now mode
extractFrames('filename')	save individual frames from a movie